

# Products for Pharmaceutical Labeling

FLX ID	Film and Topcoat	Film Type				Adhesive	Liner	Temperature				Sterilization Method								Peel Adhesion (oz/in)				Security		Printing		Quick-Ship			
		Polyester	Polypropylene	Polyethylene	Diacetate			Min Application Temperature F	Min Service Temperature F	Max Service Temperature F	Storage at 70°F - 2 Years, 50% RH	No Sterilization	Autoclave	Ethylene Oxide (ETO)	Gamma	Electron Beam (EB)	No Sterilization	Autoclave	Ethylene Oxide (ETO)	Gamma	Electron Beam (EB)	Stainless Steel	Acrylic	Glass	Polypropylene	Tamper-Evident Film	Fluorescing Adhesive		Conventional	Thermal Transfer	Digital
												Glass Vials (8mm)				HDPE Syringes (16mm)															
210	PHARMcal® PM 150 Clear TC-249	•				V-95	150 Poly	35	-40	302	•	•	•	•	•	•	•	•	•	•	42	49	42	6			•				
368	PHARMcal® PM 150 Clear TC-249	•				V-29F	150 Poly	35	-40	302	•										45	63	51	10		•	•				
212	PHARMcal® PP 200 H Clear TC-332		•			V-232	150 Poly	35	-40	176	•	•	•	•	•	•	•	•	•	•	36	44	39	6			•				
211	PHARMcal® PP 200 H Clear TC-332		•			V-232C	150 Poly	35	-40	176	•	•	•	•	•	•	•	•	•	•	36	44	39	6		•	•			•	
26253	PHARMcal® PP20041D1 - PP 200 H Clear DigiPRO®		•			V-232C	150 Poly	35	-40	176	•	•	•	•	•	•	•	•	•	•	36	44	39	6		•				•	
163	Classics Plus® 20085 - DA 200 Clear PTP				•	V-23	50# Kraft	35	-40	212	•	•									56	64	64	10	•		•			•	
214	FLEXmark® DA 200 Clear PTP				•	V-88	50# Kraft	35	-40	212	•	•									85	85	85	80	•		•			•	
27516	TAMPERmark™ PE 200 White Destruct TC-634			•		V-88	50# Kraft	35	-40	176	•	•	•	•	•	•					Tear	Tear	Tear	Tear	•		•	•		•	
443	TAMPERmark™ MM 150 Silver Void TC-249	•				L-156	50# Kraft	35	-40	302	•	•									Tear	Tear	Tear	Tear	•		•				
446	TAMPERmark™ PM 200 White Void II TC-249	•				V-344	50# Kraft	35	-40	302	•										Tear	Tear	Tear	Tear	•		•			•	
Cryogenic																															
300	CryoFLEX™ PEWM PE 350 F White Matte PTP			•		V-232	50# Kraft	35	-112	176	•	•	•	•	•	•	•	•	•	•	64	69	61	21			•	•		•	
301	CryoFLEX™ PPC - PP 200 H Clear TC-679		•			V-232	50# Kraft	35	-112	176	•	•	•	•	•	•	•	•	•	•	36	44	39	6			•	•		•	
302	CryoFLEX™ PPW - PP 230 H White TC-679		•			V-232	50# Kraft	35	-112	176	•	•	•	•	•	•	•	•	•	•	41	50	45	3			•	•		•	

These products are recommendations only and should be tested in the actual conditions of use. The typical performance data and recommendations for using FLEXcon products shall be used only as a guide. All purchasers of FLEXcon products shall be responsible for independently determining the suitability of the material for the purpose of its intended use.

250-500 MSI minimum order on most products.  
 \* Quick-Ship orders ship within 2 days.  
 5/2/2016

*because* ( Sometimes the margin for error is zero. )